

## EN • CLINI FORD

**In vitro diagnostic medical device – CE mark - In compliance with EC Directive 98/79**

Cod. AD19139P CLINI FORD 10 TEST: n.1 foil pouch containing 10 C1 cuvettes; n.1 foil pouch containing 10 S1 tubes; n.1 foil pouch containing 10 S2 tubes; n.1 GREEN capped vial containing 50µl capillaries; n.1 sachet containing disposable pipette tips; package insert.

### Intended Use

Reagent pack for the quantitative determination of the total Antioxidant Status (FORD) on whole blood, with Clini5 instruments series. For assessment of the oxido-reductive index (REDOX INDEX), when tested in combination with the FORT test (Cod. AD19107P and/or AD19109P). Clini5 instruments series is an in vitro diagnostic system intended for health care professionals.

### Reagent composition

|                                |                               |
|--------------------------------|-------------------------------|
| <i>WHITE tube - S1 Reagent</i> | <i>BLUE tube – S2 Reagent</i> |
| Hyperosmolar buffer            | pH 5.2 acetate buffer         |
| <i>Vial - S3 Reagent</i>       | <i>Cuvette - C1 Reagent</i>   |
| Iron solution                  | Chromogen                     |

### Reagent Preparation and Storage

Reagents are ready to use. Reagents are stable if stored at room temperature (15-30°C/59-86°F) and kept in the closed aluminium foil pouch up to the date marked on the packaging.

### Performance Characteristics

#### Linearity

0.80-5.00 mmol/l trolox eq.

When the reading obtained is outside the linearity range, <X or >Y is displayed, (X marks the lower end and Y the upper end).

#### Repeatability

The analytical repeatability as within-run precision was established by assaying whole blood samples and it is expressed as a percentage of the Coefficient of Variability (% CV).

| Test (n) | Mean (mmol/l trolox eq) | Std Dev | %CV  |
|----------|-------------------------|---------|------|
| 20       | 1.22                    | 0.048   | 3.89 |

#### Precision

The between series analytical precision was established by assaying blood samples and it is expressed as percent of the Coefficient of Variability (% CV).

| Test (n) | Mean (mmol/l trolox eq) | Std Dev | %CV  |
|----------|-------------------------|---------|------|
| 20       | 0.92                    | 0.039   | 4.24 |

# callegari

## Method comparison (accuracy)

A comparison study with capillary blood as specimen type measured using the Clini5 instrument instrument gave the following results:

|                           |                            |
|---------------------------|----------------------------|
| Sample number (n)         | 106                        |
| Measurement range         | 0.72-2.22 mmol/l trolox eq |
| Passing-Bablok regression | $y=1.0000x-0.0050$         |
| Correlation coefficient   | 0.934                      |
| Mean bias % (95% CI)      | -0.72 (-1.52 a +0.08)      |